

Application of wrist-worn raw accelerometry in large-scale population-based studies: challenges and field-work innovations




Kirsten Rennie, Kate Westgate, Katrien Wijndaele & Soren
Brage

NIHR Cambridge Biomedical Research Centre – Diet, Anthropometry and
Physical Activity Group

MRC Epidemiology Unit, University of Cambridge



Objective measure of physical activity in large-scale studies



Diurnal activity profile by urban-rural





Challenges and solutions are rarely in academic publications


How to translate protocols into action



Focus on adapting and implementing in LMICs





Challenges and some solutions







Infrastructure

- Customs
- Electricity supply to charge devices
- Internet access
- Deploying and retrieving devices





Training

- Large number of staff
- Limited training time allocated
- Non-specialist
- Key champions
- Standard set up of monitors
- Teaching technical quickly

Activity monitor



- Monitor worn on **non-dominant** wrist
- Wear for 7-days, 24/hours per day
- Continue with normal activities
- Waterproof
- Engraving on strap is always on left hand side of wrist





Cultural Perceptions



- Tracking or monitoring for other reason
- Identify sensitivities early on before study begins
- FAQs adapted for local settings







Fieldwork Practicalities




- Acceptability and compliance
- Wearing continuously
- Straps fitting
- Health & Safety work regulations
- Lost monitors
- Stolen monitors



Fieldwork practicalities







Data handling


- Linking device to participant
- Download of data
e.g. if collecting at 100Hz for 7 days that is 60 million rows of data to download
- Data transfer and storage





Data processing

- Ensuring high quality data
- Detecting anomalies
- Processing large datasets
- Bespoke pipeline
- Post-processing methods



Pampro Manager Current project: Fenland_1_GA Logged in as: pa.test Change project Log out

Find files Update files Run process by status Project settings Analysis settings Visualisation settings Error logs Reset process

Select	Raw data files (total = 2100)	Monitor	QC Diagnostics	Convert to HDF5	Extract Still Bouts	Calibrate Monitor	Standard Analysis
<input type="checkbox"/>	501677T_01455...	14551	Completed	Completed	Completed	Completed	Completed
<input type="checkbox"/>	507002V_left wris...	11502	Completed	Completed	Completed	Completed	Completed
<input type="checkbox"/>	507467Z_01150...	11500	Completed	Completed	Completed	Completed	Completed
<input type="checkbox"/>	509197J_left wrist...	14990	Completed	Completed	Completed	Completed	Completed
<input type="checkbox"/>	509546S_left wris...	11493	Completed	Completed	Completed	Completed	Completed
<input type="checkbox"/>	509657J_left wrist...	11461	Completed	Completed	Completed	Completed	Completed
<input type="checkbox"/>	509657V_left wris...	11459	Completed	Completed	Completed	Completed	Completed

Many thanks to:



PA Technical Team:

Lewis Griffiths, Stef Hollidge, Ella Hutchinson & Antonia Smith



Field Teams from:

University of Witwatersrand, South Africa

University of Yaounde, Cameroon

KEMRI, Kenya

Kisumu Data Collection Team

STEPS Field Team, Malawi

World Health Organisation, Geneva

Federal University of Pelotas, Pelotas, Brazil

OpenLab, Newcastle

Activinsights Ltd

<http://www.mrc-epid.cam.ac.uk/pa-seminar/>
Kirsten.rennie@mrc-epid.cam.ac.uk

